

Product datasheet for **HK200406**

Trefoil Factor 3 (TFF3) Human qPCR Template Standard (BC017859)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017859
Synonyms:	HITF; hP1.B; intestinal trefoil factor; ITF; secretory protein; TFI; trefoil factor 3; trefoil factor 3 (intestinal); trefoil factor 3, HITF, human intestinal trefoil factor
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206786) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)