

Product datasheet for **HK200412**

Transmembrane 4 L6 family member 1 (TM4SF1) Human qPCR Template Standard (BC034145)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC034145
Synonyms:	H-L6; L6; M3S1; membrane component, chromosome 3, surface marker 1; TAAL6; transmembrane 4 L six family member 1; transmembrane 4 superfamily member 1; tumor-associated antigen L6
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210707) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)