

Product datasheet for **HK200418**

TMEM16A (ANO1) Human qPCR Template Standard (BC033036)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC033036
Synonyms:	anoctamin 1, calcium activated chloride; FLJ10261; membrane protein; oral cancer overexpressed 2; ORAOV2; TAOS2; TMEM16A; transmembrane protein 16A; transmembrane protein 16A (eight membrane-spanning domains); tumor amplified and overexpressed sequence 2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212787) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

