

## Product datasheet for **HK200421**

### TMLHE Human qPCR Template Standard (BC025269)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC025269
Synonyms:	BBOX2; butyrobetaine (gamma), 2-oxoglutarate dioxygenase (gamma-butyrobetaine hydroxylase) 2; epsilon-trimethyllysine 2-oxoglutarate; epsilon-trimethyllysine hydroxylase; FLJ10727; OTTHUMP00000024255; TMLH; trimethyllysine hydroxylase, epsilon; XAP130
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212908) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)