

Product datasheet for **HK200425**

Tumor protein D52 like 3 (TPD52L3) Human qPCR Template Standard (BC017589)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017589
Synonyms:	hD55; MGC26757; MGC45374; NYD-SP25; OTTHUMP00000021042; OTTHUMP00000021043; OTTHUMP00000021044; protein kinase NYD-SP25; tumor protein D52-like 3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP233259) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)