

Product datasheet for **HK200426**

TRAK1 Human qPCR Template Standard (BC015922)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC015922
Synonyms:	OGT(O-Glc-NAc transferase)-interacting protein 106 KDa; OGT(O Glc NAc transferase) interacting protein 106 KDa; OIP106; trafficking protein, kinesin binding 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226395) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)