

## Product datasheet for **HK200427**

### TREM2 Human qPCR Template Standard (BC032362)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC032362
Synonyms:	OTTHUMP00000016368; OTTHUMP00000039771; TREM-2; Trem2a; Trem2b; Trem2c; triggering receptor expressed on monocytes 2; triggering receptor expressed on myeloid cells 2; triggering receptor expressed on myeloid cells 2a
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213304) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)