

## Product datasheet for **HK200429**

### TTC30B Human qPCR Template Standard (BC033795)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC033795
Synonyms:	2510042P03Rik; novel tetratrico peptide repeat (TPR) containing protein; OTTMUSP00000016153; Tcc30b; tetratricopeptide repeat domain 30B
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP223916) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)