

Product datasheet for **HK200431**

TTYH1 Human qPCR Template Standard (BC019358)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC019358
Synonyms:	4930459B04Rik; 6330408P11Rik; OTTMUSP00000024583; tty; tweety 1; tweety homolog 1 (Drosophila)
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225893) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)