

Product datasheet for **HK200440**

UBE2S Human qPCR Template Standard (BC066948)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC066948
Synonyms:	E2-EPF; E2EPF; EPF5; ubiquitin-conjugating enzyme E2-24 kD; ubiquitin-conjugating enzyme E2S; ubiquitin-protein ligase S; ubiquitin carrier protein S
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210937) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)