

Product datasheet for **HK200442**

UBXN8 Human qPCR Template Standard (BC020694)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC020694
Synonyms:	D8S2298E; REP8; Reproduction/chromosome 8; reproduction 8; UBXD6; UBX domain containing 6; UBX domain protein 8
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP223579) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)