

Product datasheet for **HK200451**

WDR42A (DCAF8) Human qPCR Template Standard (BC013107)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC013107
Synonyms:	DKFZp781G1096; FLJ35857; H326; MGC99640; MGC117276; MGC118891; OTTHUMP00000029883; OTTHUMP00000029884; OTTHUMP00000029885; WD repeat domain 42A
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP211805) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)