

Product datasheet for **HK200457**

YIF1B Human qPCR Template Standard (BC014974)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC014974
Synonyms:	9430029K10Rik; FinGER8; MGC130116; OTTMUSP00000030234; Yip1 interacting factor homolog B; Yip1 interacting factor homolog B (<i>S. cerevisiae</i>); Yip1b
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP232690) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)