

Product datasheet for **HK200464**

BS69 (ZMYND11) Human qPCR Template Standard (BC034784)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC034784
Synonyms:	adenovirus 5 E1A binding protein; bone morphogenetic protein receptor-ass; BRAM1; BS69; BS69 variant 1; BS69 variant 2; BS69 variant 3; BS69 variant 4; MGC111056; OTTHUMP00000018935; OTTHUMP00000018936; OTTHUMP00000018937; OTTHUMP00000044555; RP11-486H9.1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209539) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)