

Product datasheet for **HK200573**

Activin Receptor Type IA (ACVR1) Human qPCR Template Standard (NM_001105)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001105
Synonyms:	ACTRI; ACVR1A; ACVRLK2; ALK2; FOP; SKR1; TSRI
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230471) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)