

Product datasheet for **HK200645**

Angiotensin II Type 1 Receptor (AGTR1) Human qPCR Template Standard (NM_000685)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000685
Synonyms:	AG2S; AGTR1B; AT1; AT1AR; AT1B; AT1BR; AT1R; AT2R1; HAT1R
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225777) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)