

## Product datasheet for **HK200793**

### **APE1 (APEX1) Human qPCR Template Standard (NM\_080649)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>Accession No:</b>	NM_080649
<b>Synonyms:</b>	APE; APE1; APEN; APEX; APX; HAP1; REF1
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)