

Product datasheet for **HK200868**

ARMC8 Human qPCR Template Standard (NM_015396)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_015396 |
| Synonyms: | GID5; HSPC056; S863-2; VID28 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP233673) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)