

Product datasheet for **HK201010**

B3GALNT1 Human qPCR Template Standard (NM_033167)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_033167
Synonyms:	B3GALT3; beta3Gal-T3; galT3; Gb4Cer; GLCT3; GLOB; P; P1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP227418) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)