

Product datasheet for **HK201682**

Chemokine Receptor D6 (ACKR2) Human qPCR Template Standard (NM_001296)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001296
Synonyms:	CCBP2; CCR9; CCR10; CMKBR9; D6; hD6
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205211) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)