

Product datasheet for **HK201919**

Liver Carboxylesterase 1 (CES1) Human qPCR Template Standard (NM_001266)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001266
Synonyms:	ACAT; CE-1; CEH; CES2; hCE-1; HMSE; HMSE1; PCE-1; REH; SES1; TGH
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP231082) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)