

Product datasheet for **HK201960**

CHMP1B Human qPCR Template Standard (NM_020412)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_020412
Synonyms:	C10orf2; C18-ORF2; C18orf2; CHMP1.5; hVps46-2; Vps46-2; Vps46B
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213669) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)