

Product datasheet for **HK201968**

CHN 1 (CHN1) Human qPCR Template Standard (NM_001822)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001822
Synonyms:	ARHGAP2; CHN; DURS2; NC; RHOGAP2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP229991) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)