

## Product datasheet for **HK202011**

### Creatine kinase B type (CKB) Human qPCR Template Standard (NM\_001823)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001823
Synonyms:	B-CK; BCK; CKBB; CPK-B; HEL-211; HEL-S-29
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205604) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)