

Product datasheet for **HK202236**

MCSF Receptor (CSF1R) Human qPCR Template Standard (NM_005211)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_005211
Synonyms:	BANDDOS; C-FMS; CD115; CSF-1R; CSFR; FIM2; FMS; HDLS; M-CSF-R
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208389) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)