

## Product datasheet for **HK202364**

### Aromatase (CYP19A1) Human qPCR Template Standard (NM\_031226)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_031226
Synonyms:	ARO; ARO1; CPV1; CYAR; CYP19; CYPXIX; P-450AROM
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200089) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)