

Product datasheet for **HK202366**

Cytochrome P450 1A2 (CYP1A2) Human qPCR Template Standard (NM_000761)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000761
Synonyms:	CP12; CYPIA2; P3-450; P450(PA)
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200705) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)