

Product datasheet for **HK202370**

Cytochrome P450 2B6 (CYP2B6) Human qPCR Template Standard (NM_000767)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000767
Synonyms:	CPB6; CYP2B; CYP2B7; CYP2B7P; CYP11B6; EFVM; IIB1; P450
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200710) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)