

## Product datasheet for **HK202373**

### Cytochrome P450 2E1 (CYP2E1) Human qPCR Template Standard (NM\_000773)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000773
Synonyms:	CPE1; CYP2E; P450-J; P450C2E
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200715) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)