

Product datasheet for **HK202376**

Cytochrome P450 3A4 (CYP3A4) Human qPCR Template Standard (NM_017460)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_017460
Synonyms:	CP33; CP34; CYP3A; CYP3A3; CYP11A3; CYP11A4; HLP; NF-25; P450C3; P450PCN1; VDDR3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212375) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)