

## Product datasheet for **HK202522**

### **DHRS9 Human qPCR Template Standard (NM\_199204)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_199204
<b>Synonyms:</b>	3-alpha-HSD; 3ALPHA-HSD; RDH-E2; RDH-TBE; RDH15; RDHL; RDHTBE; RETSDR8; SDR9C4
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP208848) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)