

Product datasheet for **HK202532**

DATF1 (DIDO1) Human qPCR Template Standard (NM_080797)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_080797
Synonyms:	BYE1; C20orf158; DATF-1; DATF1; DIDO2; DIDO3; DIO-1; DIO1; dj885L7.8
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP234707) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)