

## Product datasheet for **HK202548**

### DLEU2 Human qPCR Template Standard (NM\_006021)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_006021
Synonyms:	1810047A16Rik; 4932411A06Rik; 7730401J12Rik; AI197277; AI788824; Alt1; deleted in lymphocytic leukemia, 2; DLB2, LEU2, BCMSUN, MGC88287; Leu2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209045) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)