

## Product datasheet for **HK202602**

### DPAGT1 Human qPCR Template Standard (NM\_001382)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001382
Synonyms:	ALG7; CDG-lj; CDG1j; CMS13; CMSTA2; D11S366; DGPT; DPAGT; DPAGT2; G1PT; GPT; UAGT; UGAT
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205274) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)