

Product datasheet for **HK202704**

EEF1A1 Human qPCR Template Standard (NM_001402)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001402
Synonyms:	CCS-3; CCS3; EE1A1; EEF-1; EEF1A; eEF1A-1; EF-Tu; EF1A; GRAF-1EF; LENG7; PTI1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205289) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)