

## Product datasheet for **HK202757**

### EIF4A1 Human qPCR Template Standard (NM\_001416)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001416
Synonyms:	DDX2A; EIF-4A; eIF-4A-I; EIF4A; eIF4A-I
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205300) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)