

## Product datasheet for **HK202759**

### EIF4A2 Human qPCR Template Standard (NM\_001967)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001967
Synonyms:	BM-010; DDX2B; eIF-4A-II; EIF4A; eIF4A-II; EIF4F
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205728) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)