

## Product datasheet for **HK202846**

### GTP binding protein era homolog (ERAL1) Human qPCR Template Standard (NM\_005702)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_005702
Synonyms:	CEGA; ERA; ERA-W; ERAL1A; ERAL1B; H-ERA; HERA-A; HERA-B; PRLTS6
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208793) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)