

Product datasheet for **HK202850**

XPG (ERCC5) Human qPCR Template Standard (NM_000123)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000123
Synonyms:	COFS3; ERCC5-201; ERCM2; UVDR; XPG; XPGC
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200108) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)