

## Product datasheet for **HK202852**

### ERGIC3 Human qPCR Template Standard (NM\_015966)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_015966
Synonyms:	C2orf47; C20orf47; CGI-54; dj477O4.2; Erv46; NY-BR-84; PRO0989; SDBCAG84
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225421) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)