

Product datasheet for **HK202903**

Junctional Adhesion Molecule 1 (F11R) Human qPCR Template Standard (NM_144504)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_144504
Synonyms:	JAM, KAT, JAM1, JAMA, JCAM, CD321, JAM-1, JAM-A, PAM-1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP212339) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)