

Product datasheet for **HK203213**

FUSIP1 (SRSF10) Human qPCR Template Standard (NM_054016)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_054016
Synonyms:	FUSIP1; FUSIP2; NSSR; PPP1R149; SFRS13; SFRS13A; SRp38; SRp40; TASR; TASR1; TASR2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226809) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)