

Product datasheet for **HK203529**

Grainyhead like protein 1 homolog (GRHL1) Human qPCR Template Standard (NM_014552)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_014552
Synonyms:	grainyhead-like 1 (Drosophila); LBP-32; LBP32; LBP protein 32; leader-binding protein 32; mammalian grainyhead; MGR; TFPC2L2; transcription factor CP2-like 2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210962) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)