

Product datasheet for **HK203607**

MACROH2A1 Human qPCR Template Standard (NM_138610)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_138610
Synonyms:	H2A.y; H2A/y; H2AF12M; H2AFY; MACROH2A1.1; macroH2A1.2; mH2A1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226986) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)