

## Product datasheet for **HK203654**

### HDAC1 Human qPCR Template Standard (NM\_004964)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004964
Synonyms:	GON-10; HD1; KDAC1; RPD3; RPD3L1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208195) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)