

## Product datasheet for **HK203690**

### HIF-1 alpha (HIF1A) Human qPCR Template Standard (NM\_001530)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001530
Synonyms:	bHLHe78; HIF-1-alpha; HIF-1A; HIF-1alpha; HIF1; HIF1-ALPHA; MOP1; PASD8
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205393) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)