

Product datasheet for **HK203706**

H2BC8 Human qPCR Template Standard (NM_003518)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003518
Synonyms:	dj221C16.8; H2B.1A; H2B/a; H2BC4; H2BC6; H2BC7; H2BC10; H2BFA; HIST1H2BG
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207022) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)