

Product datasheet for **HK203710**

H2AC18 Human qPCR Template Standard (NM_003516)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003516
Synonyms:	H2A; H2a-615; H2A.2; H2A/O; H2A/q; H2AC19; H2AFO; HIST2H2AA; HIST2H2AA3
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207020) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)