

Product datasheet for **HK203711**

H2BC21 Human qPCR Template Standard (NM_003528)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003528
Synonyms:	GL105; H2B; H2B.1; H2BE; H2BFQ; H2BGL105; H2BQ; HIST2H2BE
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207032) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)