

## Product datasheet for **HK203713**

### Hexokinase 1 (HK1) Human qPCR Template Standard (NM\_000188)

#### Product data:

|                  |   |
|------------------|---|
| Product Type:    | qPCR Template Standards   |
| ACCN:            | NM_000188   |
| Synonyms:        | hexokinase; HK; HK1-ta; HK1-tb; HK1-tc; HKD; HKI; HMSNR; HXK1; NEDVIBA; RP79  |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (HP229819) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| Storage:         | Store at -20°C.   |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.   |



[View online »](#)